

Gostim2 User Guide

Christophe Pallier

Contents

| | |
|--|----------|
| Gostim2 User Guide | 2 |
| What is Gostim2? | 2 |
| GUI vs. CLI: which should I use? | 2 |
| Installation | 3 |
| Pre-built binaries (recommended) | 3 |
| Building from source | 3 |
| The input file: TSV/CSV format | 3 |
| Required columns | 3 |
| Extra columns | 4 |
| Minimal example | 4 |
| Important notes on timing | 4 |
| Asset directories | 4 |
| Stimulus types | 5 |
| IMAGE | 5 |
| SOUND | 5 |
| TEXT | 5 |
| BOX | 5 |
| IMAGE_STREAM | 6 |
| TEXT_STREAM | 6 |
| SOUND_STREAM | 6 |
| VIDEO | 7 |
| Running the GUI (<code>gostim2-gui</code>) | 7 |
| Running the CLI (<code>gostim2</code>) | 7 |
| Required argument | 7 |
| File and directory options | 8 |
| Display options | 8 |
| Font options | 8 |
| Fixation cross options | 9 |
| Timing options | 9 |
| Other options | 9 |
| Example invocations | 9 |
| Output logs | 10 |
| Log structure | 10 |
| Results directory resolution | 10 |
| Advanced topics | 10 |
| Linux: maximum timing precision | 10 |

| | |
|--|----|
| Hardware triggers (DLP-IO8-G) | 11 |
| Config file (TOML) | 11 |
| Quick reference: stimulus type summary | 11 |
| License and credits | 12 |

Gostim2 User Guide

Version: 0.2.0 **Author:** Christophe Pallier **License:** GNU General Public License v3

What is Gostim2?

Gostim2 is a multimedia stimulus delivery system designed for experimental psychology and cognitive neuroscience. It presents images, sounds, text, and video according to a **fixed schedule defined in a plain-text table** — no programming required.

Gostim2 is best suited for experiments where:

- Stimuli are presented in a predetermined order and at known times.
- Precise timing is critical (fMRI, MEG, EEG studies).
- You want a simple, portable setup that runs on Linux, Windows, and macOS.

Important limitation: Gostim2 does **not** support real-time branching or feedback. The presentation schedule is fully determined before the experiment starts. All key-press events are recorded, but they cannot alter stimulus delivery. If you need adaptive or interactive paradigms, consider [goxpyriment](#) (Go) or [expyriment](#) (Python).

GUI vs. CLI: which should I use?

Gostim2 ships as **two separate programs** that share the same engine:

| | <code>gostim2-gui</code> (graphical) | <code>gostim2</code> (command-line) |
|-----------------------------|---|---|
| How to start | Double-click or run from terminal | Run from a terminal with flags |
| Configuration | Interactive form — point-and-click | Flags and/or a TOML config file |
| Settings persistence | Saved automatically to <code>gostim2_config.toml</code> | Must be re-specified each time (or use <code>-config</code>) |
| Best for | First-time users, occasional experimenters | Scripting, automation, HPC clusters, TTY-only servers |
| Requires display | Yes | No (can run headless in TTY mode on Linux) |
| Timing | Identical to CLI | Identical to GUI |

Both programs use **the same timing engine** and produce **identical output logs**. Choosing one over the other is purely a matter of workflow preference.

Installation

Pre-built binaries (recommended)

Download the appropriate package for your platform from the [GitHub Releases page](#):

| Platform | Package |
|---------------------|---|
| Windows | <code>gostim2-setup.exe</code> (installer) or <code>.zip</code> |
| macOS Apple Silicon | <code>.app.zip</code> or <code>macos-arm64.tar.gz</code> |
| macOS Intel | <code>macos-x86_64.tar.gz</code> |
| Linux x86_64 | <code>.AppImage</code> or <code>.deb</code> or <code>.tar.gz</code> |

macOS note: The first time you open the app, macOS may show a security warning because the binary is not signed. See [macOS installation and security](#) for step-by-step instructions to bypass this.

Building from source

Requires Go 1.25+ and the SDL3 libraries (`sd13`, `sd13_image`, `sd13_ttf`):

```
git clone https://github.com/chrplr/gostim2
cd gostim2
make build           # builds gostim2 and gostim2-gui for the current platform
make install        # copies binaries to $GOPATH/bin
```

The input file: TSV/CSV format

The entire experiment is described in a single **tab-separated values (TSV)** or **comma-separated values (CSV)** file. You can create and edit this file in any spreadsheet program (LibreOffice Calc, Excel) or a plain-text editor.

The delimiter is detected automatically (tab, comma, semicolon, or pipe).

Required columns

Your file **must** contain these four column headers (case-insensitive):

| Column | Type | Description |
|-------------------------|--------------|--|
| <code>onset_time</code> | integer (ms) | When the stimulus should appear, in milliseconds from experiment start |
| <code>duration</code> | integer (ms) | How long to display/play the stimulus (see per-type notes below) |

| Column | Type | Description |
|----------------------|--------|---|
| <code>type</code> | string | Stimulus type — one of the values described in the next section |
| <code>stimuli</code> | string | The content: a filename, a text string, or a stream specification |

Extra columns

Any additional columns (e.g. `condition`, `trial_id`, `expected_response`) are ignored during playback but are **preserved verbatim in the output log**, making it easy to merge logs with your design matrix.

Minimal example

```
onset_time duration type stimuli
1000 500 IMAGE face01.png
2000 300 TEXT Hello, world!
4000 2000 BOX Please press\nany key
7000 1 SOUND beep.wav
```

Important notes on timing

- All times are in **milliseconds**.
- `onset_time` is measured from the moment the participant presses the key to start the experiment.
- Stimuli must not overlap in time. The validator will reject a file with overlapping onsets.
- For `SOUND` stimuli, set `duration` to 1 (or any positive integer) — the sound plays to its natural end regardless of this value.

Asset directories

By default, Gostim2 looks for media files (images, sounds) in a folder named `stimuli/` or `assets/` located **next to the CSV file**. You can override this with the `-stimuli-dir / -assets` flag (CLI) or the corresponding field in the GUI.

A recommended project layout:

```
my_experiment/
  run1.csv
  stimuli/
    face01.png
    face02.png
    beep.wav
  gostim2-results/ (created automatically)
```

Stimulus types

IMAGE

Displays a single image file for the specified duration.

```
type      stimuli
IMAGE     face01.png
```

Supported formats: PNG, JPEG, BMP, GIF (first frame only), and other formats supported by `SDL3_image`.

SOUND

Plays a single audio file. The sound plays to its **natural end** regardless of the `duration` value; set `duration` to 1 as a placeholder.

```
onset_time  duration  type      stimuli
7000        1         SOUND     beep.wav
```

Supported formats: **WAV**, **FLAC**, **MP3**, **OGG**. Audio is automatically resampled if the sample rate does not match the device output rate.

TEXT

Renders a single line of text centered on the screen.

```
type      stimuli
TEXT     Fixation cross coming up...
```

Font, size, text colour, and background colour are all configurable (see [Display options](#)).

BOX

Renders **multiline text** centered on the screen. Use the two-character sequence `\n` within the `stimuli` string to insert a line break.

```
type      stimuli
BOX       Please press\nthe SPACE bar\nto continue
```

Renders as:

```
Please press
the SPACE bar
to continue
```

IMAGE_STREAM

Plays a rapid sequence of images, one after another, within a single row of the CSV. Images are separated by ~ in the `stimuli` column.

Simple form — all frames use the same duration from the `duration` column, with no inter-frame gap:

| onset_time | duration | type | stimuli |
|------------|----------|--------------|----------------------------------|
| 2000 | 200 | IMAGE_STREAM | face01.png~face02.png~face03.png |

Timed form — each frame specifies its own duration and gap (blank screen between frames), using the format `filename:duration:gap`:

| onset_time | duration | type | stimuli |
|------------|----------|--------------|--|
| 2000 | 1 | IMAGE_STREAM | face01.png:200:50~face02.png:200:50~face03.png:200:0 |

Here each face is shown for 200 ms followed by a 50 ms blank screen, except the last which has no trailing gap.

TEXT_STREAM

Rapid sequence of text strings. Uses the same ~ separator and optional `:duration:gap` timing as IMAGE_STREAM.

| onset_time | duration | type | stimuli |
|------------|----------|-------------|------------------------|
| 5000 | 300 | TEXT_STREAM | The~cat~sat~on~the~mat |

Each word is displayed for 300 ms (from the `duration` column) with no gap.

With explicit timing:

| onset_time | duration | type | stimuli |
|------------|----------|-------------|----------------------------------|
| 5000 | 1 | TEXT_STREAM | The:300:50~cat:300:50~sat:300:50 |

SOUND_STREAM

Plays a sequence of audio files. Uses the same ~ separator and `:duration:gap` format. For audio streams, `duration` acts as the **SOA (Stimulus Onset Asynchrony)** — the time between the start of one sound and the start of the next.

| onset_time | duration | type | stimuli |
|------------|----------|--------------|-------------------------------|
| 3000 | 500 | SOUND_STREAM | word1.wav~word2.wav~word3.wav |

Each sound starts 500 ms after the previous one (regardless of the sound's actual length).

With explicit per-item timing:

| onset_time | duration | type | stimuli |
|------------|----------|--------------|---|
| 3000 | 1 | SOUND_STREAM | word1.wav:500:0~word2.wav:700:0~word3.wav:500:0 |

VIDEO

Plays a video file. If `duration` is set to 1, the video plays to its natural end; set a value in ms to cut it off early.

| onset_time | duration | type | stimuli |
|------------|----------|-------|----------|
| 10000 | 1 | VIDEO | clip.mp4 |

Running the GUI (`gostim2-gui`)

Launch `gostim2-gui` (double-click the icon, or run it from the terminal). A setup window appears:

1. **Experiment CSV** — click “...” to browse for your `.csv` or `.tsv` file.
2. **Stimuli directory** — usually auto-detected as the `stimuli/` or `assets/` folder next to the CSV; override if needed.
3. **Subject ID** — enter a participant identifier (always starts empty, never pre-filled).
4. **Results directory** — where output logs are saved (default: `gostim2-results/` next to the CSV).
5. **Resolution** — pick a fixed resolution or check **Autodetect resolution** to use your monitor’s native resolution.
6. **Window mode:**
 - *Windowed*: standard window.
 - *Fullscreen Desktop*: borderless window at desktop resolution; OS compositor remains active. Safest option.
 - *Fullscreen Exclusive*: takes exclusive control of the display, bypasses the compositor. Lowest latency; recommended for EEG/MEG/fMRI.
7. **Font / font size, colours** (background, text, fixation cross).
8. **VSync** — enabled by default; disable only for VRR/G-Sync setups.
9. Click the green **START** button.

When the experiment window opens, press any key to begin (unless **Skip wait** is checked). Press **Escape** at any time to abort.

Settings are automatically saved to `gostim2_config.toml` in the same directory, and reloaded at the next launch.

Running the CLI (`gostim2`)

```
gostim2 -csv experiment.csv [options]
```

Required argument

| Flag | Description |
|--------------------------------|---|
| <code>-csv <path></code> | Path to the stimulus CSV/TSV file. <code>-tsv</code> is an alias. |

File and directory options

| Flag | Default | Description |
|---|------------------------------|---|
| <code>-stimuli-dir <path></code> | auto-detected | Directory containing image/sound assets. <code>-assets</code> is an alias. |
| <code>-results-dir <path></code> | <code>gostim2-results</code> | Directory for output log files (resolved relative to the CSV file). |
| <code>-subject <id></code> | unknown | Participant identifier, included in the output filename and log. |
| <code>-config <path></code> | — | Load parameters from a TOML config file (same format as <code>gostim2_config.toml</code>). |
| <code>-start-splash <path></code> | — | Image to display before the experiment (waits for key press). |
| <code>-end-splash <path></code> | — | Image to display after the experiment (waits for key press). |

Display options

| Flag | Default | Description |
|--|-----------------|---|
| <code>-res <WxH \ Autodetect></code> | — | Screen resolution, e.g. 1920x1080. Use <code>Autodetect</code> for native resolution. |
| <code>-width <int></code> | 1440 | Screen width (ignored if <code>-res</code> is set). |
| <code>-height <int></code> | 900 | Screen height (ignored if <code>-res</code> is set). |
| <code>-display <int></code> | 0 | Index of the monitor to use. |
| <code>-scale <float></code> | 1.0 | Logical scale factor applied to all stimuli. |
| <code>-fullscreen</code> | off | Exclusive fullscreen (bypasses compositor). |
| <code>-fullscreen-desktop</code> | off | Borderless fullscreen (compositor active). |
| <code>-bg-color <R,G,B,A></code> | 0,0,0,255 | Background colour as RGBA values 0–255. |
| <code>-text-color <R,G,B,A></code> | 255,255,255,255 | Text colour as RGBA values 0–255. |
| <code>-fixation-color <R,G,B,A></code> | 255,255,255,255 | Fixation cross colour. |

Font options

| Flag | Default | Description |
|-------------------------------------|----------------------|--------------------------|
| <code>-font <path></code> | built-in Inconsolata | Path to a TTF font file. |
| <code>-font-size <int></code> | 50 | Font size in points. |

Fixation cross options

| Flag | Default | Description |
|-------------------------------|---------|--|
| <i>(default)</i> | | Fixation cross shown only on blank screens. |
| <code>-no-fixation</code> | | Never display the fixation cross. |
| <code>-fixation-always</code> | | Always display the fixation cross (superimposed on stimuli). |

Timing options

| Flag | Default | Description |
|------------------------|---------|--|
| <code>-vsync</code> | on | Synchronize rendering to the display refresh rate. |
| <code>-no-vsync</code> | | Disable VSync. |
| <code>-vrr</code> | off | Variable Refresh Rate mode: disables VSync and uses busy-wait for ~2 ms precision. For G-Sync/FreeSync monitors. |

Other options

| Flag | Description |
|----------------------------------|---|
| <code>-skip-wait</code> | Skip the “Press any key to start” prompt. |
| <code>-dlp <device></code> | Serial device path for DLP-IO8-G hardware triggers (e.g. <code>/dev/ttyUSB0</code>). |
| <code>-version</code> | Print version information and exit. |

Example invocations

Basic run

```
gostim2 -csv study1/run1.csv -subject sub-01
```

Fullscreen exclusive at 1920x1080, custom font

```
gostim2 -csv run1.csv -fullscreen -res 1920x1080 -font ~/fonts/DejaVuSans.ttf -subject 42
```

```
# Load most settings from a config file, override subject only
gostim2 -config my_experiment.toml -subject sub-07
```

```
# Run from a TTY console (Linux) for minimal latency
gostim2 -csv run1.csv -fullscreen -subject sub-01
```

Output logs

After each run, Gostim2 writes a timestamped TSV file to the results directory. The filename follows the pattern:

```
{csv_basename}_sub-{subject_id}_{YYYY-MM-DDTHH-MM-SS}.tsv
```

Log structure

The file starts with a **metadata header** (lines beginning with #) recording system information: SDL version, audio/video drivers, platform, display mode, and the exact flags used.

This is followed by a data table with these columns:

| Column | Description |
|----------------------|--|
| subject_id | The participant identifier from <code>-subject</code> or the GUI field. |
| intended_ms | The onset time as specified in the CSV. |
| actual_ms | The real onset time measured at runtime. |
| type | Event type: <code>IMAGE_ONSET</code> , <code>SOUND_ONSET</code> , <code>TEXT_ONSET</code> , etc., or <code>RESPONSE</code> . |
| label | The stimulus filename or text string. |
| <i>extra columns</i> | All extra columns from the original CSV are reproduced here. |

Results directory resolution

The results path is always resolved **relative to the CSV file**, not to the working directory:

| Command | Results saved in |
|---|--------------------------------------|
| <code>gostim2 -csv study1/run1.csv</code> | <code>study1/gostim2-results/</code> |
| <code>gostim2 -csv study1/run1.csv -results-dir data</code> | <code>study1/data/</code> |
| <code>gostim2 -csv study1/run1.csv -results-dir ""</code> | <code>study1/</code> |

Advanced topics

Linux: maximum timing precision

On Linux, you can achieve the lowest possible latency by running Gostim2 from a bare TTY console (Ctrl+Alt+F3), after stopping the display manager:

```
sudo systemctl stop gdm # or lightdm, sddm, etc.
gostim2 -csv run1.csv -fullscreen -subject sub-01
```

This allows SDL3 to use the **Direct Rendering Manager (DRM)** directly, bypassing Wayland/X11 entirely.

Hardware triggers (DLP-IO8-G)

For EEG/MEG setups requiring TTL triggers, connect a DLP-IO8-G USB device and pass its serial port:

```
gostim2 -csv run1.csv -dlp /dev/ttyUSB0 -subject sub-01
```

On Windows the port is typically COM3 or similar.

Config file (TOML)

The GUI automatically writes `gostim2_config.toml` after each session. You can hand-edit this file or create one from scratch and pass it with `-config`:

```
CsvFile       = "run1.csv"
StimuliDir    = "stimuli"
ResultsDir    = "gostim2-results"
FullscreenMode = "exclusive"
Resolution    = "1920x1080"
FontSize      = 60
BgColor       = [0, 0, 0, 255]
TextColor     = [255, 255, 255, 255]
```

Quick reference: stimulus type summary

| Type | Content of <code>stimuli</code> column | duration |
|--------------|---|--|
| IMAGE | Single image filename | Display time (ms) |
| SOUND | Single audio filename | Use 1 (plays to end) |
| TEXT | Single line of text | Display time (ms) |
| BOX | Multiline text; <code>\n</code> for line breaks | Display time (ms) |
| IMAGE_STREAM | <code>img1.png~img2.png</code> or <code>img1.png:dur:gap~...</code> | Default frame duration if no per-item timing |
| TEXT_STREAM | <code>word1~word2</code> or <code>word1:dur:gap~...</code> | Default frame duration if no per-item timing |
| SOUND_STREAM | <code>snd1.wav~snd2.wav</code> or <code>snd1.wav:soa:gap~...</code> | Default SOA if no per-item timing |
| VIDEO | Single video filename | 1 to play to end, or cut-off time (ms) |

License and credits

Gostim2 is a port of [audiovis](#) to Go, using [go-sdl3](#) bindings.

Author: [Christophe Pallier](#) christophe@pallier.org

Source code: <https://github.com/chrplr/gostim2>

Distributed under the [GNU General Public License v3](#).